

VIDAWELL

# HEATING PAD

---

## OPERATION AND STORAGE MANUAL

Read in full and retain for future use



Keep this manual with the appliance. If the appliance is to be used by a third party, this instruction manual must be supplied with it.

**MODEL: VWHP-01**

# HEATING PAD

## FEATURES

---

- 3 temperature settings
- Machine washable
- Soothing heat therapy
- Runtime of 90 minutes
- Auto cutoff
- Fast heat up
- Overheat protection
- Fleece material

You can use the heat pad to gently warm and relax various parts of the body. Using the heat pad stimulates blood circulation in hardened muscles and helps you to relax after a strenuous day. The heat pad features an electric temperature control that regulates the temperature according to the chosen setting.

## BEFORE FIRST USE

---

- Before first use, remove any stickers or similar from the appliance except for the rating label.
- Dispose of any plastic bags as they pose a risk of suffocation to young children.

Please read this manual before operating the appliance and keep it in a safe place for future reference. If this appliance is transferred to a new owner, the manual should also be provided to the new owner.

## IMPORTANT SAFETY MEASURES AND PRECAUTIONS

---

- Carefully check the heat pad before each use, if you notice wear or damage.
- Do **not** use if you notice wear, damage or signs of improper use on the heat pad, the switch or the cables. Return it to the supplier.
- Before connecting the device to your power supply, please ensure that the supply voltage stated on the rating plate is compatible with your mains supply.
- Do not fold or crease the heat pad during use.
- Neither pierce the pad with, nor attach, safety pins or other spiky or sharp objects to the pad.
- Do not use the heat pad on children or persons who are disabled, asleep or sensitive to heat.
- Children should be supervised to ensure that they do not play with the appliance.

- Never use the pad unsupervised.
- Do not fall asleep while the heat pad is switched on.
- Prolonged use of the heated cushion at a high setting could lead to skin burns.
- Do not cover the pad with any other cushions. The control switch should not be covered, or placed on or below the pad while the unit is in operation.
- Do not sit on the heat pad. Instead, rest the pad on or against the part of the body to be treated.
- Never touch a heat pad that has fallen into water. Unplug the unit from the mains outlet immediately.
- Keep the mains lead away from hot surfaces.
- Never carry, pull, or turn the heat pad by the mains lead and do not let the lead get tangled.
- Do not use the heat pad when wet, and only use it in dry environments (not in the bathroom or similar environments).
- The controls and leads must not be exposed to any kind of moisture.
- If a fault occurs, do not attempt to repair the unit yourself. Repairs must only be carried out by an authorized specialist dealer or other suitably qualified personnel.
- A damaged power supply cable must be replaced with a power cable of the same type.
- When storing the heat pad, allow it to cool down before folding.
- Should you have any concerns about health issues, consult your doctor before using the heat pad.
- Do not use the heat pad on parts of the body that are swollen, injured or inflamed.
- Should you experience persistent pains in the muscles or joints, please inform your doctor. Persistent pains can be a symptom of serious illness.
- If you experience any pain or discomfort while using the unit, stop using it immediately.

## **WARNINGS**

- Please ensure that the polythene packing is kept away from the reach of children. Risk of suffocation.
- Do not puncture the heated cushion.
- Do not use the heated cushion when it is folded.
- Only use the heated cushion indoors.
- The heating pad can be washed by hand. You may wash the cushion cover in a normal wash at max. 30 °C.
- Do not chlorinate.
- Do not dry the heated cushion in a tumble dryer. The cushion cover may not be put in a dryer.
- Do not iron the heated cushion. You may iron the cushion cover at a low temperature.
- Do not dry clean.

## OPERATING INSTRUCTIONS

---

1. Plug the unit into the mains outlet and move the slider control from position 0 to position 1.
2. The power LED will light up in green. After a few minutes you will be able to feel the pad getting noticeably warmer. If you require a higher temperature setting, move the switch to the next position, 2, or to the maximum heat setting, position 3. If the pad gets too hot for you, move the slider back to position 2 or 1.
3. After a short while the heating pad cools down noticeably.
4. To switch the unit off, move the switch back to position 0.
5. The green power LED will go out, indicating that the unit is switched off. Unplug the unit from the mains outlet. Move the selector switch to position 1 or 2 before falling asleep or any continuous use. The heat pad will switch off automatically after approx. 90 minutes of continuous use.
6. To switch it back on again, move the slider control to position O and then back to the desired heat setting.
7. Unplug the unit from the mains outlet once you have finished using it.

## CLEANING AND MAINTANANCE

---

- Before cleaning the pad, unplug it from the mains outlet and let it cool down for at least ten minutes.
- The heating pad is equipped with a detachable cable. Disconnect the plug and remove the cable from the heating pad.
- You may wash the heating pad gently by hand. The best way to do this is to put the heating pad in a bath with lukewarm water and some mild detergent and then gently squeeze it out.
- Rinse the heating pad several times to remove all the detergent.
- Store the unit in a clean and dry place, laid out flat with nothing covering or resting on top of it.
- Only use the appliance once it is completely dry.

## DISPOSAL

---

This product must not be disposed of together with domestic waste.

## POWER RATING

---

220-240V-50HZ, 100W

**Customer Support**

[www.asiapacbrands.com](http://www.asiapacbrands.com)

Helpline New Zealand call toll free 0508 274 2722

Australia call +64 9 974 9206

**Distributed by**

Asia Pacific Brands 2019 Ltd

PO Box 301145, Albany, Auckland 0752, New Zealand

**[www.vidawell.net](http://www.vidawell.net)**